




# **SECURE CHAIN OF CUSTODY PROCEDURES FOR THOROUGHBRED RACING TEST SPECIMENS**

**BEFORE VIEWING THIS PRESENTATION, THE VISTOR IS  
ENCOURAGED TO CLICK THE RED BUTTON AT THE BOTTOM  
OF THE “ABOUT US” SECTION ON THE HOME PAGE**

# WHAT MAKES A SIGNAKEY SPECIAL

- Highly encrypted 2-D code – not in the public domain.
- Every Signakey is unambiguously unique.
- Easy to decode with IOS or Droid based Smartphones.
- App is free - but must be registered to each device.
- User must be approved by Commission.
- Highly secure, permissions based, fixed length blockchain.



# **KHRC AND LESSON LEARNED**

**Dr. Bruce Howard, Equine Medical Director, KHRC**

“KHRC has never won litigation for use of forbidden substances from a big purse race”

Because of:

- Inadequate chain-of-custody procedure and practices.
- Underfunded legal team - versus owners like Amr Zedan - Medina Spirit

# **BETTER C-O-C FOR TXRC AND TVMDL**


## **REQUIREMENTS:**

- KEEP ALL CURRENT PROCEDURES AND SUPPLIES.
- ADD THREE NEW TAMPER EVIDENT SEAL BAGS PER HORSE
- VETS WILL NEED IOS DEVICES WITH SIGNAKEY RACE-TRACK APP
- LAB AND SPC WILL NEED WEBCAMS & SIGNAKEY SOFTWARE
- SIGNAKEY WILL SUPPLY APP, WEBCAMS, SOFTWARE & TRAINING
  
- PILOT DURATION 12 WEEKS

# TAMPER-EVIDENT SEAL SPECIMEN BAGS

SAMPLE ID LABEL

TAMPER EVIDENT SEAL



BAG #1


UNIQUE SIGNAKEY

BAG #2


BAG #3

# STEPS FOR SCANNING A SIGNAKEY


## Open APP - GoTo Settings for Profile




## Enter Credentials Enable Phone GPS



## Center Object and App will Scan SignaKey



## Results Summary after Scanning SignaKey




# IN THE TEST BARN

- USE ALL EXISTING PROCEDURES TO COLLECT 1 URINE & THEN 4 BLOOD SAMPLES.
- PLACE ALL SAMPLES IN TAMPER EVIDENT BAG #1 AND SECURE SEAL.
- USING iPhone WITH RACE-TRACK APP – SCAN SIGNAKEY ALL WHILE IN THE STALL AREA.
- **THE SCAN BEGINS THE SECURE CHAIN-OF-CUSTODY FOR THE SPECIMENS.**
- ATTACH THE ID STICKER TO BAG #1 AND MOVE IT TO THE SAMPLE PREPARATION CENTER.

# TE BAG #1 WITH LARGE URINE SAMPLE & 4 VIALS

TAMPER EVIDENT SEAL BAG WITH PREPRINTED TEXT COVERED BY LABELS FOR DEMONSTRATION




SAMPLE # LABEL

SIGNAKEY LABEL




# IN THE SAMPLE PREPARATION CENTER

- CHECK BAG #1 SEAL IS SECURE. GET PHOTO.
- IF YES – SCAN SIGNAKEY USING COMPUTER WEB-CAM
- **NEXT STEP IN THE CHAIN-OF-CUSTODY PROCESS.**
- OPEN BAG #1 REMOVE CONTENTS.
- CENTRIFUGE VIALS - SPLIT URINE IN 2 EQUAL SAMPLES
- FOLLOW ALL ESTABLISHED SAMPLE PREP PROCEDURES.
- 3 BLOOD & 1 URINE SAMPLES GO INTO THE LAB BAG
- 1 BLOOD & 1 URINE SAMPLE GO INTO SPLIT BAG
- SECURE TE SEAL ON EACH BAG - SCAN BOTH SIGNAKEY.
- **NEXT STEP IN THE CHAIN OF CUSTODY PROCESS**
- PUT LAB BAG IN SHIPPING CONTAINER IN LOCKING FRIDGE
- SPLIT SAMPLE BAGS GOES INTO SECURE LOCKING FREEZER.




# IN SPC - FILL IDENTICAL SAMPLE BAGS



GET **4** BLOOD VIALS



GET **2** URINE SAMPLES




**TAKE POST-RACE SPECIMENS  
IN TEST BARN AND PLACE IN  
THE TAMPER EVIDENT BAG #1**



**FROM BARN**



**LOGIN WITH SMARTPHONE TO  
ESTABLISH THE USER**



**SCANNING AT EACH LOCATION  
CAPTURES TIMESTAMP &  
GEOSPATIAL COORDINATES**

**TO LAB**

**CONFIRM SEAL INTACT THEN  
SCAN SIGNAKEY BEFORE  
OPENING BAG**



# UN-CUSTOMIZED DATABASE TEMPLATE

File Action Help Thursday, 24 Mar 2022 15:15:56 UTC-04:00 Configure TXRC.Bill Log Out

Network  Camera  Decode  Record  Camera Web Tries Database: Texas Racing Commision Test Tracking Element: Case

Read **Create** Update Delete Verify SignalKey History

Key	<input checked="" type="radio"/> New: <input type="radio"/> Decode <b>Idle</b> Clear
Status	<input checked="" type="radio"/> New: <input type="radio"/>
Venue	<input checked="" type="radio"/> New: <input type="radio"/>
Geolocation	<input checked="" type="radio"/> New: <input type="radio"/>
<b>Horse</b>	
Name	<input checked="" type="radio"/> New: <input type="radio"/>
Image	<input checked="" type="radio"/> New: <input type="radio"/> Upload Browse Browse to File
Registration	<input checked="" type="radio"/> New: <input type="radio"/>
Image	<input type="radio"/> New: <input type="radio"/> Upload Browse Browse to File
<b>Primary</b>	
Seal	<input type="radio"/> New: <input type="radio"/> Decode <b>Idle</b> Clear
Lab	<input type="radio"/> New: <input type="radio"/>
Geolocation	<input type="radio"/> New: <input type="radio"/>
Report	<input type="radio"/> New: <input type="radio"/> Upload Browse Browse to File
<b>Split</b>	
Seal	<input type="radio"/> New: <input type="radio"/> Decode <b>Idle</b> Clear
Lab	<input type="radio"/> New: <input type="radio"/>
Geolocation	<input type="radio"/> New: <input type="radio"/>
Report	<input type="radio"/> New: <input type="radio"/> Upload Browse Browse to File

Reset Do Create Cancel

08:16:43 SignaKey Cloud Client, Version 1.12.41.0, © 2022 SignaKey  
08:16:44 Initialization complete.  
08:17:09 Manually logged in as TXRC.Bill  
08:17:14 Opened database Texas Racing Commision Test Tracking, owned by TXRC

BAG #1 ID

TRACK DATA

HORSE DATA

PRIMARY BAG #2

SPLIT BAG #3

# AT THE LAB

- ✓ CHECK THAT BAG SEAL IS SECURE – TAKE CONFIRMING PHOTO.
- ✓ IF NO – REFUSE DELIVERY.
- ✓ IF YES – **SCAN SIGNAKEY WITH WEB-CAM – CONFIRMS C-O-C**
- ✓ OPEN SAMPLE BAG – REMOVE SPECIMENS – SAVE SAMPLE BAG.
- ✓ GENERATE ACCESSION NUMBER AND CREATE SAMPLE STICKERS
- ✓ PERFORM ALL TESTS USING NORMAL PRESCRIBED METHODS.
- ✓ OPEN RACE-TRACK SOFTWARE
- ✓ SCAN SIGNAKEY ON SAMPLE BAG WITH WEBCAM.

# COMPLETE RECORD FORMAT

File Action Help Friday, 25 Mar 2022 08:42:25 UTC-04:00 Configure TXRC.Bill Log Out

Network Camera Decode Record Camera Tries Database: Element:  
Texas Racing Commision Test Tracking Case

Read Create Update Delete Verify SignaKey History

+ Add a filter or read all Case elements.  Counts Only Fields

SignalKey Cloud Client Read Results

Export 1 of Case

Key	Status	Venue	Geolocation	Horse				Primary			Split				
				Name	Image	Registration	Image	Seal	Lab	Geolocation	Report	Seal	Lab	Geolocation	Report
1-130-1-1	Failed	Houston Race	29.76045, -95.36978	Trigger		12345		1-130-2-1000	TVMDL			1-130-2-1001	Split Test Lab		

SCAN OF SIGNAKEY

UPLOAD RESULTS

Reset Do Read Cancel

08:16:43 SignalKey Cloud Client, Version 1.12.41.0, © 2022 SignalKey  
08:16:44 Initialization complete.  
08:17:09 Manually logged in as TXRC.Bill  
08:17:14 Opened database Texas Racing Commision Test Tracking, owned by TXRC  
08:40:23 • Read Success • Found 1 match.

# AT THE LAB

- ✓ SCAN OF LAB BAG SIGNAKEY BRINGS UP CORRECT FILE.
- ✓ SCAN ALSO INFORMS SOFTWARE OF PENDING STATUS CHANGE.
- ✓ UPLOAD ALL TEST RESULTS FROM LIM INTO RACE-TRACK.
- ✓ IF ALL TESTS ARE NEGATIVE – DISCARD SAMPLE BAG.
- ✓ IF TESTS ARE POSITIVE – RETAIN BAG WITH CASE FILE EVIDENCE.
- ✓ SELECTED SPLIT LAB TO BE INFORMED BY TXRC OF REQUIREMENTS.

# IN OOC SETTING OPTION

- ✓ SCHEDULE WITH TRAINER & TXRC VET OR VET TECH.
- ✓ NEED TE SAMPLE BAG/4 VIALS/URINE CUP/HAIR COMB & BAG.
- ✓ CHECK FOR IPHONE SIGNAL COVERAGE.
- ✓ USE ALL SAME PROCEDURES AS BARN PLUS HAIR IF NO URINE.
- ✓ PUT SAMPLES IN BAG, SECURE SEAL, SCAN SIGNAKEY.
- ✓ COMPLETE SIGNATURE CARD WITH TRAINER WITNESS.
- ✓ RETURN BAG TO SPC, SCAN SAME AS RACE DAY FOR C-O-C.



# COMPLETE RECORD FORMAT

File Action Help Friday, 25 Mar 2022 08:42:25 UTC-04:00 [Configure](#) [TXRC.Bill](#) [Log Out](#)

Network 
Camera 
Decode 
Record 
Camera Web
Tries
Database: Texas Racing Commission Test Tracking
Element: Case

Read Create Update Delete Verify SignaKey History

+ Add a filter or read all Case elements.  Counts Only Fields

◆ SignaKey Cloud Client Read Results

[Export](#) 1 of Case

				Horse				Primary				Split			
Key	Status	Venue	Geolocation	Name	Image	Registration	Image	Seal	Lab	Geolocation	Report	Seal	Lab	Geolocation	Report
1-130-1-1	Failed	Houston Race	29.76045, -95.36978	Trigger		12345		1-130-2-1000	TVMDL		<a href="#">PDF</a>	1-130-2-1001	Split Test Lab		<a href="#">PDF</a>



**DRUG TESTING LABORATORY**  
480 Agronomy Road, TAMU College Station, TX, (866) 646-5623

## CASE FILE

Customer	TX Racing Commission	Accession #	C202160001
Sample ID	RP045528	Date Collected	07/30/2020
Matrix	Urine	Date Received	08/03/2020
Drug(s) Confirmed	<b>Morphine</b>	Date Reported	08/25/2020

Jim Blodgett  
 Director of Investigations  
 Texas Racing Commission  
 8505 Cross Park Drive, Suite 110  
 Austin, TX 78754

Received: 2020-09-03  
 Laboratory No.: EACL-200903-4  
 External ID: RP045528  
 Sample No: **411136**  
 No. of Samples: 1

INVESTIGATION: Split Sample: Morphine  
 One urine sample submitted from the Texas Racing Commission, RP045528, was delivered on September 3<sup>rd</sup> 2020.

The urine sample was subjected to qualitative analysis for the presence of morphine and morphine-3-glucuronide. Detection and identification were performed utilizing liquid chromatography - mass spectrometry. Morphine and morphine glucuronide were detected in the sample. No further testing has been assigned.

If you have any questions or require additional information, please don't hesitate to contact me.

The remainder of the original sample will be stored at the Maddy Lab and disposed of after 2 months.

# STATUS SUMMARY

- SAMPLE BAGS ARE GENERIC – SHOULD BE TXRC CUSTOMIZED.
- DATA BASE WORKS – BUT NEEDS TO BE CUSTOMIZED.
- PHONE APP “RACE-TRACK” DESIGNED BUT NOT YET FULLY PROGRAMED
- 12 WEEK PILOT THRU JUNE 59 RACE DAYS - APR 22 THRU JUNE 30.
- EVERYTHING DESIGNED TO BE HISA COMPLIANT